



# LiquiBand® Exceed™

## Drawing the line on even, consistent wound coverage

LiquiBand® Exceed™ a **2-octyl cyanoacrylate adhesive** has a unique felt, applicator tip to facilitate even, consistent application.




### ✓ MARKET LEADING YIELD

- ▶ 0.8ml per applicator allows for closure of wounds up to 30cm in length
- ▶ No clog technology allows intraoperative reuse for up to 90 minutes on a single patient
- ▶ Single applicator can be used for most closures

### ✓ PURE OCTYL FORMULATION

- ▶ Maximises durability and flexibility
- ▶ Provides high viscosity
- ▶ Creates intimate skin contact



### ✓ UNIQUE APPLICATOR TIP

- ▶ Elliptical shape allows for narrow or wide application
- ▶ Porous felt tip allows drip-free priming and even, consistent application

### ✓ EASY TO USE APPLICATOR

- ▶ Winged applicator allows for safe and easy activation
- ▶ Provides controlled delivery of adhesive





Better wound closure starts with a better tool




## LiquiBand® Exceed™

Offers strong and secure wound closure equivalent to market leading products



## Key Benefits

- ✔ **Better Yield of Adhesive<sup>1</sup>**
  - ▶ Sufficient adhesive for up to 30cm of wound coverage
  - ▶ Two times the yield of market leading competitor
- ✔ **Unique No Clog Tip<sup>1</sup>**
  - ▶ 0.8mL volume provides industry leading coverage with a single applicator
  - ▶ Suitable for intraoperative reuse for up to 90 min on a single patient
- ✔ **Controlled Expression of Adhesive<sup>1</sup>**
  - ▶ Applicator offers discrete control during expression of adhesive
- ✔ **Microbial Barrier<sup>1</sup>**
  - ▶ Acts as a barrier against gram-positive, gram-negative and fungi microbes. Proven effectiveness against S. aureus, P.aeruginosa, E.coli, Candida albicans and MRSA
- ✔ **Maintains wound closure throughout healing process<sup>1</sup>**
  - ▶ The adhesive will stay on the skin for 5-10 days and slough off naturally



**T** 1-855-GLUES-77  
**E** info@liquiband.com  
**W** liquiband.com